

Double Side- Chemical Immersion Dunk Tank.

BSL3/BSL4 Pathogen Neutralising Solution



- Lockable access covers with compression seal
- Wall plate containment bulkhead for facility interface
- Large 750 x 500 vision panel (APR standard)
- 316 Stainless construction
- Complete with submersion paddles and drain plug



PBSC Bio decontamination chambers are in use throughout the world – Years of expertise have been distilled down into a cost effective Dunk Tank package. Easy to install, easy to validate

For more information contact:

PBSC Ltd. Bradley Junction Industrial Estate – Leeds Road – Huddersfield UK

Tel +44 (0)1484 354500 Email info@pbsc.co.uk Web www.pbsc.co.uk

